

## Inertsil ODS-2

- Virtually free of trace metal contaminants which can cause peak tailing
- Outstanding peak shape for acids and bases
- Exceptional physical and chemical durability provides long column life
- Vast database of chromatographic applications available to assist in method development
- One of most popular and trusted HPLC columns in industry for nearly 15 years

The original 2-Series Inertsil phases, including ODS-2, C8, Ph, C4, and SIL-150A°, created the high purity silica revolution in the late 1980s. These phases utilize a 99.999% pure silica gel which is extremely durable and free from dents and cracks which can cause premature column failure.

The original Inertsil 2-Series of columns remain some of the most useful and popular on the market today. However, GL Sciences' Inertsil 3-Series phases were a major advancement on the original 2-Series, and generally provide superior chromatography and alternative selectivity to the 2-Series. We recommend starting out with Inertsil ODS-3 or Inertsil ODS-3V for new methods using an ODS (C18) phase.

Particle Size(s)	Particle Shape	Surface Area	Pore Size	Pore Volume	Silica Purity	Bonded Phase	End-Capping	Carbon Load	pH Range*
5 µm	Spherical	320 m <sup>2</sup> /g	150 Å	1.15 mL/g	99.999%	Octadecyl Groups	Yes	18.5%	2 - 7.5

[www.winlab.com.au](http://www.winlab.com.au)

Ph: 61 07 3205 5233

Fax: 61 07 3205 1209