


WINLAB PTY LTD

2 Pinnacle St, Brendale, Queensland, Australia, 4500
PO Box 5007, Brendale, Queensland, Australia, 4500

Ph: +61 7 3205 5233: Fax: +61 7 3205 1209

Email: info@winlab.com.au www.winlab.com.au



Histology – Cytology – Hangers

Hangers are made from Polyacetal and Stainless Steel, durable and resistant to solvents, having proved themselves over many years in laboratories around the world.



AW101-3-1 Replacement part, compatible with all SR30 Slide Racks on Linear and X-Y staining machines.

AW101-3-2 Replacement part, compatible with all SR20 Slide Racks on Linear and X-Y staining machines.

AW101-3-1M Replacement part, compatible with Instrumec 30 place Slide Racks on Linear and X-Y staining machines.

AW101-3-2M Replacement part, compatible with Instrumec 19 place Slide Racks on Linear and X-Y staining machines.

AW101-3-9 Dual Stainer Double Carrier compatible with Instrumec 19 Place Racks

AW101-3-17 Dual Stainer Single Carrier compatible with Instrumec 19 Place Racks



AW101-3-10 Dual Stainer Double Carrier compatible with Winlab SR20 Racks

AW101-3-18 Dual Stainer Single Carrier compatible with Winlab SR20 Racks

AW101-3-11 Dual Stainer Carrier compatible with Winlab SR30 Racks



AW101-3-14 SR20 Slide Rack Hanger compatible with Medite Stainers

AW101-3-16 Instrumec 19 Slide Rack Hanger compatible with Medite Stainers

AW101-3-12 SR20 Series Plastic Handle, economical, microwave safe, and solvent resistant for manual staining, to avoid chemical & heat contact.

AW101-3-13 SR30 Series Plastic Handle, economical, microwave safe, and solvent resistant for manual staining to avoid chemical & heat contact.

