


**WINLAB PTY LTD**

2 Pinnacle St, Brendale, Queensland, Australia, 4500

PO Box 5007, Brendale, Queensland, Australia, 4500

Ph: +61 7 3205 5233: Fax: +61 7 3205 1209

Email: [info@winlab.com.au](mailto:info@winlab.com.au) [www.winlab.com.au](http://www.winlab.com.au)



## Histology and Cytology Slide Racks

**Winlab's Machine Slide Racks are made from a PolyCeramic process developed in our workshops, are extremely rigid, and highly resistant to solvents up to 200°C.**

**"Just drop one on concrete and hear the difference"**



AW101-2-1 SRLS30 A good all-rounder used by machine manufacturers; Vee Grooves for the hanger are vertical, horizontal and 45°; compatible with Leica, Tharmac, R A Lamb, and Thermo Stainers.  
Max Temps 200°C



AW101-2-2 SRM30 Exclusively made to suit the Medite X-Y Series of Stainers, having a thinner body to fit narrow troughs, and includes the swivel feature found on that type of machine.  
Max Temps 80°C



AW101-2-3 SRLSM20 For machine integration, having a narrow bottom half and a larger top half, this rack is ideally suited to the narrow troughs in the Medite X-Y Series and Sakura type Glass Coverslipping machines. Vee Grooves for the hanger are vertical, horizontal and 45°.  
Max Temps 200°C



AW101-2-4 SRM20 Exclusively made to suit the Medite X-Y Series of Stainers, having a slim body to fit narrow troughs and includes the swivel feature found on that type of machine.  
Max Temps 200°C.



AW101-2-5 SRLS20 Compatible with all Stainers and Coverslippers. Vee Grooves for the hanger are vertical, horizontal and 45°.  
Max Temps 200°C



AW101-2-8 SRLM30 A good all-rounder used by machine manufacturers, having a swivel action hanger compatible with Medite, Leica, Tharmac, R A Lamb, and Thermo Stainers.  
Max Temps 200°C